

PTA & INJECTION. COMBINED.

Chameleon™ PTA Balloon Catheter



SERVES MULTIPLE PURPOSES

Perform an entire procedure with both angioplasty and injection of diagnostic or therapeutic fluids



REDUCES NEED FOR CATHETER EXCHANGE

Perform angiograms directly through the Chameleon™ balloon catheter



ENABLES HANDS-FREE REFLUX ANGIOGRAPHY

Perform reflux angiography while the Chameleon™ balloon catheter occludes venous outflow, which eliminates manual compression and enhances operator safety



HELPS OPTIMIZE USE OF CONTRAST

Precise delivery of contrast directly to the treatment site may reduce the amount of contrast needed compared to standard PTA




MINIMIZES RADIATION EXPOSURE

Fewer device exchanges and increased operator distance from radiation source may result in lower radiation exposure¹

The Chameleon™ PTA balloon catheter uniquely combines the functionality of a high pressure balloon catheter and a diagnostic catheter — which enables proximal injection of diagnostic fluids (like contrast media) and therapeutic fluids.

Versatility and efficiency

The Chameleon™ PTA balloon catheter enables you to visualize while injecting through the catheter — whether the balloon is inflated or deflated — all while maintaining wire position.




Promising Clinical Results with the Chameleon™ PTA Balloon Catheter


CASE 1

Courtesy of Dr. Jeffrey Hoggard, Chapel Hill, North Carolina
Juxta-anastomotic stenosis of AV fistula
Chameleon™ PTA Balloon Catheter 8 x 40mm


RETROGRADE APPROACH




Juxta-anastomotic narrowing



Balloon inflation



Balloon inflation



Final imaging through the side port of the Chameleon™ balloon catheter


CASE 2

Courtesy of Dr. Robert Madden, West Springfield, Massachusetts
Upper Arm Loop Graft, Antegrade
Chameleon™ PTA Balloon Catheter 8 x 40mm


REFLUX ANGIOGRAPHY




Balloon inflation



Balloon inflation



Hands-free reflux angioplasty through the inflated Chameleon™ balloon catheter



Outflow imaging through the deflated Chameleon™ balloon catheter

Available in a range of sizes to cover most dialysis access procedures.

Catalog Number	Diameter (mm)	Balloon Length (mm)	Catheter Length (cm)	Sheath Size (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Guidewire (in)
CH05-40-75US	5	40	75	6	14	25	0.035
CH06-40-75US	6	40	75	6	14	25	0.035
CH07-40-75US	7	40	75	6	14	22	0.035
CH08-40-75US	8	40	75	6	14	22	0.035
CH09-40-75US	9	40	75	7	12	18	0.035
CH10-40-75US	10	40	75	7	12	14	0.035
CH12-40-75US	12	40	75	8	12	14	0.035

1. Crawford J, Kokkosis A, Kim P, Gasparis A, Labropoulos N. Fistuloplasty using a radiation-and-time-saving sheathless balloon catheter. J Vasc Access. 2019;20(3):276-280.

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 11/2019 – US180606(1) – [WF# 3623282]