




Argyle™

safety fistula cannula with anti-reflux valve



1. PUNCTURE NEEDLE
2. SIDE HOLES
3. CANNULA
4. CONE
5. LUER LOCK FITTING
6. SAFETY CAP
7. TAB
8. BEVEL MARK
9. NEEDLE HUB
10. V-PLUG

PRODUCT INFORMATION

Item code	Cannula			Color
	Size I.D. (mm)	Size O.D. (mm)	Length (mm)	
1538VSVS	1.53	1.8 (15 G)	38	Blue-Gray
1530VSVS	1.53	1.8 (15 G)	30	Blue-Gray
1638VSVS	1.3	1.6 (16 G)	38	White
1630VSVS	1.3	1.6 (16 G)	30	White
1738VSVS	1.12	1.4 (17 G)	38	Red-Violet
1730VSVS	1.12	1.4 (17 G)	30	Red-Violet

FLOW TABLES

Length mm (in)	15G		16G		17G	
	Flow Rate (mL/min)	Arterial Pressure (mmHg)	Flow Rate (mL/min)	Arterial Pressure (mmHg)	Flow Rate (mL/min)	Arterial Pressure (mmHg)
38 mm (1.50 in)*	450	-214	300	-221	200	-188

Arterial pressure should not exceed care provider guidelines. Results shown above achieved using 3-4 cP Glycerin / Water solution. The extracorporeal circuit includes devices in addition to the fistula cannula; observe the lowest flow rate for all devices within the system. * The longest cannula was tested to represent the worst case scenario.

MATERIAL INFORMATION

Needle	Stainless steel
Cannula	Polypropylene
Cone	Polypropylene

SPECIAL FEATURES

- Sharp metal needle integrating a safety mechanism designed to prevent needle-stick injuries
- Anti-reflux valve to prevent blood backflow
- Blunt flexible cannula with 4 side holes
- Cannula hole information:
15 and 16 gauge, all styles have hole size of 0.65 mm, with tolerances of (0.2 mm, -0.35 mm)
17 gauge, all styles, have hole size of 0.50 mm with tolerances of (0.1 and -0.2 mm)

BIOCOMPATIBILITY

Latex-free	✓
DEHP-free	✓

PACKAGING

	Units	Material
Product packaging (pouch)	1	The packaging film is a polyethylene / nylon blend
Department packaging (box)	50	Paper cardboard
Transport packaging (carton)	300	Paper cardboard

SHELF LIFE AND STERILIZATION

Sterilization method	E-Beam radiation
Expiration time	2 years

Medtronic

DISCLAIMER: Specifications and other information presented in this document are for information purpose only and do not provide any recommendation or advise of any kind. Data may differ from actual values quoted here and Medtronic does not make any warranties or representations to any characteristics. Subject to changes without notice.

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.
™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
EMEA-RC-2000017 - 20-weu-argyle-fistula-cannula-technical-sheet-4742126

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)