



**EMODIAL**

**Medtronic**

exclusive distributor in Canada

**EMOSTAN<sup>®</sup>**

Fistula Pad



**EMODIAL**

Via Finati, 4/A - Ferrara (ITALY)  
Tel. (+39) 0532 978198 - Fax (+39) 0532 978414  
[www.emodial.it](http://www.emodial.it) - [info@emodial.it](mailto:info@emodial.it)



REV. 11/2018



**Description:**


Highly absorbent pad.

**Indications:**

Compressive haemostatic action, specific for post-dialysis arterovenous access care to be used after the removal of the fistula needles. Also suitable for patients who require intravenous treatment and blood transfusions.

**Purpose:**

Designed to stop bleeding and to reduce the haemorrhagic risk thanks to the compression on the needles holes and to work as a barrier against the risk of infection. The non-adherent, biocompatible layer on the pad doesn't stick to the skin **avoiding the clot removal and the risk of bleeding.**



## Fistula Pad

**Directions for use**

- 1** Put the pad over the needle site and hold in place with thumb.



- 2** Remove the needle and press for a suitable time to stop bleeding.



REFERENCES Data and technical specifications can be changed without notice

CODE	PAD DIMENSION cm	PAD THICKNESS cm	Pcs/ Pack	Packs/ Box
S0210-1	ø 2.6	0.8 ± 0.1	1	700
S0211-1	ø 3.5	0.8 ± 0.1	1	700
S0212-1	3,5 x 4.5	0.8 ± 0.1	1	500

**Classification:**  
**Biocompatible:**  
**Sterilization:**  
**Expiration:**



Class I medical device  
according to UNI EN ISO 10993 standards  
Gamma rays  
5 years  
Single-use